


# Guide to Non-Coring Needles

All implanted ports must be accessed with a non-coring needle, sometimes referred to as a Huber Needle. Non-coring needles are manufactured by many companies with each having specific features. You will need to have knowledge of the brand used in your institution and reference the manufacturer’s instructions for use. Non-coring needles come in many sizes to fit the depth of the port and viscosity of the infusate. Several companies have developed power injector ports used for injecting contrast materials in enhanced imaging procedures. These power ports must be accessed with a power non-coring and extension. Always use safety engineered devices for needlestick prevention. Below is an overview of the safety non-coring needle devices available on the market.

[Kawasumi Laboratories America, Inc Product Catalog](#)

Photo	Product
	<p><a href="#">K Shield® Power Injectable</a> <a href="#">Video</a></p>
	<p><a href="#">K Shield®</a></p>
	<p><a href="#">K Shield Advantage®</a></p>


## [Angiodynamics](#)

Photo	Product
	<p><a href="#">LifeGuard Safety Huber Needle</a> <a href="#">LifeGuard Promotional Literature PDF</a></p>






**These lists are for reference and educational use only. Pedagogy does not recommend the use of any product or product manufacturer.**

**Learn more at <http://www.pedagogyeducation.com/infusion>**

Osiris Medical

Photo	Product
	<p><a href="#">Empress Safety Huber Needle</a> <a href="#">Training Demonstration and Video</a> <a href="#">Empress/Osiris Instruction sheet PDF</a></p>

Bard Access Systems - Bard Safety Port Access Needles

Photo	Product
	<p><a href="#">EZ Huber™ Safety Infusion Set</a></p>
	<p><a href="#">HuberPlus™ Safety Infusion Set</a></p>
	<p><a href="#">LiftLoc™ Safety Infusion Set</a></p>
	<p><a href="#">MiniLoc™ Safety Infusion Set</a></p>
	<p><a href="#">SafeStep™ Huber Needle Set</a></p>


**These lists are for reference and educational use only. Pedagogy does not recommend the use of any product or product manufacturer.**

**Learn more at <http://www.pedagogyeducation.com/infusion>**

Bard Access Systems - Bard Power Injectable Port Access Needles

Photo	Product
	<a href="#"><u>Powerloc® EZ Safety Infusion Set</u></a>
	<a href="#"><u>Powerloc® Max Safety Infusion Set</u></a>
	<a href="#"><u>Powerloc® Safety Infusion Set</u></a>

B-Braun Safety Huber Needles

Photo	Product
	<a href="#"><u>B. Braun's WHIN® Safe Huber Needle</u></a>

Smiths Medical

Photo	Product
	<a href="#"><u>Gripper Plus Safety Huber Needle</u></a> <a href="#"><u>Gripper Plus Huber Instruction Sheet</u></a>
	<a href="#"><u>Gripper Micro</u></a>

**These lists are for reference and educational use only. Pedagogy does not recommend the use of any product or product manufacturer.**

**Learn more at <http://www.pedagogyeducation.com/infusion>**